




Sequências de DNA na identificação de espécies e análise filogenética de microrganismos

Dra. Chirlei Glienke

Dra. Desirrê Petters-Vandresen


Que genes usar?

- O que usar para procariotos?
 - Para gêneros: gene 16S do rRNA
 - Para espécies: Sequenciamento parcial de genes house-keeping:
 - Exemplos: genes *atpD* (ATP synthase F1, beta subunit), *gyrB* (DNA gyrase B subunit), *rpoB* (RNA polymerase beta subunit), *recA* (recombinase A) and *trpB* (tryptophan synthetase, beta subunit).



Que genes usar?


- O que usar para eucariotos?
- Para fungos:
 - Para gêneros: LSU (28S do rRNA)
 - ITS1-5.8S-ITS2 do rDNA, Segmento $\alpha 1$ do Fator de Elongação da Tradução; Sequencias parciais dos genes: β -tubulina, GPDH, Histona, Actina, rpb2 e Mating Type
- Para trypanosomas:
- glucose-6-phosphate isomerase (GPI), 18S ribosomal RNA (SSU rRNA), glycosomal GAPDH (gGAPDH), DNA Mitochondrial




SECUENCIAS DE DNA PARA FILOGENIA DE LEVEDURAS - YeastIP

Marker	Description	Position in the gene*	Length (bp)
LSU	Complete sequence of the large subunit 26S ribosomal RNA gene	65–3364	3300
D1/D2 LSU	Partial sequence of the 26S ribosomal gene comprising the D1/D2 region	65–636	570
SSU	Complete sequence of the small subunit 18S ribosomal RNA gene	1–1800	1700
ITS	Ribosomal RNA region containing the intergenic region 1 (between 18S and 5.8S), the 5.8S ribosomal RNA gene and the intergenic region 2 (between 5.8S and 26S)	First base of ITS1 to last base of ITS2	400–600
mtSSU	Mitochondrial small subunit 15S ribosomal RNA gene	383–1006	400–600
<i>RPB1</i>	Partial sequence of the RNA polymerase II largest subunit coding gene	253–873	620
<i>RBP2</i>	Partial sequence of the RNA polymerase II second largest subunit coding gene	1645–2319	680–1000
<i>TEF1</i> -alpha	Partial sequence of the translation elongation factor 1-alpha coding gene	64–1190	930
<i>ACT1</i>	Partial sequence of the exon2 of the actin coding gene	405–1383	980†
mtCOX II	Partial sequence of mitochondrial cytochrome C oxidase subunit 2 coding gene	121–707	590


(A)





DNA BARCODING

A TOOL FOR SPECIMEN IDENTIFICATION AND SPECIES
DISCOVERY



<https://ibol.org/about/dna-barcoding/>

DNA BARCODING

Como funciona?



Step 1: Isolate DNA from the sample



Step 2: Amplify the target DNA barcode region using PCR

Step 3: Sequence the PCR products



Step 4: Compare the resulting sequences against reference databases to find the matching species

<https://ibol.org/about/dna-barcoding/>



Evolution



Ecology



Collections, data, & informatics



Conservation



THE DNA BARCODE

[HOME](#)[ABOUT ▾](#)[PROGRAMS ▾](#)[NEWS & MEDIA ▾](#)[JOIN](#)


Plant barcoding studies use one or a few plastid regions (e.g. *rbcL* and *matK*, and the non-coding spacer *trnH-psbA*) and the internal transcribed spacer (ITS) region of nuclear ribosomal DNA.




Animal barcoding studies use a region in the mitochondrial cytochrome *c* oxidase 1 gene (“CO1”).



Fungal barcoding studies use the internal transcribed spacer (ITS) region in the nuclear ribosomal cistron. This region shows reasonable discriminatory power at the species level in many groups.



DNA BARCODING



THE LIBRARY

Barcode sequences are placed in the Barcode of Life Data Systems (BOLD) database – an online workbench that includes a reference library of DNA barcodes that can be used to assign identities to sequences of unknown origin.



BOLD is a searchable repository for barcode records, storing specimen data and images as well as sequences and trace files. It provides an identification engine based on the current barcode library and monitors the number of barcode sequence records and species coverage.



DNA BARCODING

Fungos: nuclear internal transcribed spacer (ITS) do RNA ribossomal –
Barcoding primário

Fusarium e Trichoderma necessidade de translation elongation factor 1
alpha (TEF-1 α) como barcode secundário

Arbuscular mycorrhizal (AM) and in rust fungi: small subunit (nuSSU)
and the large subunit (nuLSU) do RNA ribossomal RNA


Oomycota: mitochondrial cytochrome oxidase c subunits (COX1 and
COX2).

BARCODE OF LIFE DATA SYSTEM v4

**Advancing biodiversity science
through DNA-based species
identification.**

<https://www.boldsystems.org/>

QBOL – DNA Barcoding para pragas quarentenárias vegetais



QBOL

Development of a new diagnostic tool using DNA barcoding to identify quarantine organisms in support of plant health

Continue to

- About QBOL
- e-Newsletter
- Training courses
- Meeting 2012

DNA Barcoding **Activities** **Q-lists** **Databases** **Collaboration** **Partners**

QBOL project has come to an end

On 21 September 2012 the QBOL project ended officially. Many results were obtained in the different work packages.

➤ [Read more](#)

Centralizado na Europa

<https://www.qbol.org/en/qbol.htm>



[DNA Barcoding](#)

[Activities](#)

[Q-lists](#)

[Databases](#)

[C](#)

Barcoding of
Arthropods

Barcoding of Bacteria

Barcoding of Fungi

Barcoding of
Nematodes

Barcoding of
Phytoplasmata

Barcoding of Viruses

DNA Banks

Library/Database/Informatics

Validation/Evaluation

Dissemination



Q-BOL Activities

QBOL Activities:

DNA Barcoding de fungos – CBS na Holanda



DNA Barcoding **Activities** Q-lists Databases Collaboration Partners

- Barcoding of Arthropods
- Barcoding of Bacteria
- Barcoding of Fungi <**
- Barcoding of Nematodes
- Barcoding of Phytoplasmas
- Barcoding of Viruses
- DNA Banks
- Library/Database/Information
- Validation/Evaluation
- Dissemination



Barcoding of Fungi


In this work package a short list of 19 Q-species were selected for barcoding.

Contact

Ewald Groenewald
[View profile >>](#)

WP Coordinator: KNAW-CBS
(Partner 9)

KNAW-CBS, PRI, FERA and CIP

[DNA Barcoding](#)[Activities](#)[Q-lists](#)[Databases](#)[Collaboration](#)[Partners](#)[Barcoding of
Arthropods](#)[Barcoding of Bacteria](#)[Barcoding of Fungi](#)[Barcoding of
Nematodes <](#)[Barcoding of
Phytoplasmas](#)[Barcoding of Viruses](#)[DNA Banks](#)[Library/Database/Informa](#)[Validation/Evaluation](#)[Dissemination](#)

Barcoding of Nematodes

In this work package a base list of 32 nematode species was created for which barcodes needed to be collected

Contact


Juerg Frey

[View profile >>](#)


WP Coordinator: ACW (Partner 5)

ACW, PRI, VLAGIEW (ILVO).
INRA and CIP


RNA RIBOSSOMAL EUCARIÓTICO



REGIÃO ITS1-5.8S-ITS2 DO rRNA




ACTIN (1325 bp total; 1128 coding for 376 amino acids)



Exon
Intron

Primer available
New primer
Published primer

BETA-TUBULIN (1933 bp total; 1344 coding for 448 amino acids)



Exon


Intron

Primer available

New primer


Published primer

CALMODULIN (970 bp total; 450 coding for 150 amino acids)




ELONGATION FACTOR 1-ALPHA

(1975 bp total; 1383 coding for 461 amino acids)



HISTONE H3 (582 bp total; 411 coding for 137 amino acids)




Exon

Intron

Primer available
New primer
Published primer

ITS1, ITS2 & 5.8S rRNA GENE (475 bp, 158 coding)




 Coding

 Insertion

Primer available
New primer
Published primer


LARGE SUBUNIT (28S) rRNA GENE (3378 bp)



 Coding
 Insertion

Primer available
New primer
Published primer

SMALL SUBUNIT (18S) rRNA GENE (2880 bp total, 1803 bp coding)



 Coding

 Insertion

Primer available
New primer
Published primer

**Então, qual gene
vou usar?**

**Depende da
pergunta/objetivo/
problema a
resolver**






Fig. 1. Phylogenetic relationships of the Botryosphaeriales using Bayesian reconstruction and six gene portions (LSU, SSU, ITS, EF1, BT and mtSSU). Numbers above branches indicate bootstrap values/posterior probabilities. Numbers highlighted in red below branches indicate estimated dates in million years with the 95 % Highest Posterior Density interval given in brackets. Clades 1–4 in the Botryosphaeriales are indicated by a circled number on the corresponding node.





COMO FUNCIONA O SEQUENCIAMENTO?





FIGURA 1. Deoxiribonucleotídeos (dNTP) e dideoxiribonucleotídeos trifosfatados (ddNTP).


SEQUENCIAMENTO



SEQUENCIAMENTO




SEQUENCIAMENTO




SEQUENCIAMENTO

4a




Análise dos produtos fluorescentes de sequenciamento em canais separados do gel de poliacrilamida.

4b




Análise dos produtos fluorescentes de sequenciamento em um único canal do gel de poliacrilamida.

4c



Eletroforegrama do gel

SEQUENCIAMENTO



SEQUENCIAMENTO


T

C

G

```
1  gtgcccagtg agcgaggact gca
61  ggggacccga acaacatgga tag
121 acgttccatg gggccaaggga ggt
181 atgggcctca tatacaacag gat
241 tgtgccactt gtgaacagat tgc
301 accaccaatc cactaatcag gca
361 gctatggaac agatggctgg atc
421 cagactaggc agatggtaca tgc
481 ggtctgaaag atgaccttct tga
541 atgcagcgat tcaagtgatc ctc
```


select all 100 sequences selected							GenBank	Graphics	Distance tree of results	
	Description	Max Score	Total Score	Query Cover	E value	Per. Ident	Accession			
✓	Fungal endophyte isolate 522.18S ribosomal RNA gene, partial	1062	1062	100%	0.0	99.49%	K901532.1			menziesii
✓	Trametes cf. cubensis 11B-25 18S ribosomal RNA gene, partial	1062	1062	100%	0.0	99.49%	K864514.1			menziesii
✓	Basidiomycota sp. FPF30a internal transcribed spacer 1, partial	1062	1062	100%	0.0	99.49%	JX416570.1			menziesii
✓	Basidiomycota sp. FPF38a internal transcribed spacer 1, partial	1062	1062	100%	0.0	99.49%	JX416577.1			menziesii
✓	Trametes cf. cubensis 11B-20 18S ribosomal RNA gene, partial	1057	1057	100%	0.0	99.32%	KJ654515.1			menziesii
✓	Trametes cf. cubensis 8-3-5-3-B-34NFB D9-A 18S ribosomal RNA	1046	1046	100%	0.0	99.97%	KJ654513.1			menziesii
✓	Uncultured fungus clone L3042493-122-057-F05 internal transcri	1046	1046	100%	0.0	99.97%	GQ999293.1			menziesii
✓	Uncultured fungus clone L342032-122-062-107 internal transcri	1046	1046	100%	0.0	99.97%	GQ999222.1			menziesii
✓	Uncultured fungus clone L3042493-122-057-F03 internal transcri	1040	1040	100%	0.0	99.90%	GQ999291.1			menziesii
✓	Leucophaea elegans isolate PCT22 18S ribosomal RNA gene, part	1026	1026	96%	0.0	99.47%	HQ248217.1			menziesii
✓	Leucophaea lactinea voucher NMMA 2015-172 Mushrooms Ciber	1018	1018	100%	0.0	99.13%	MH010528.1			cubensis
✓	Trametes cf. cubensis 11F-21b 18S ribosomal RNA gene, partial	1018	1018	100%	0.0	99.13%	KJ654419.1			cubensis
✓	Trametes cf. cubensis 11F-21A 18S ribosomal RNA gene, partial	1018	1018	100%	0.0	99.13%	KJ654416.1			cubensis
✓	Trametes cf. cubensis 4-1-5-2-A-55(FBA)A 18S ribosomal RNA g	1016	1016	96%	0.0	99.12%	KJ654510.1			menziesii
✓	Leucophaea lactinea voucher S.D. Russell MycoAmp 6095 sma	1014	1014	100%	0.0	97.99%	MK584504.1			cubensis
✓	Trametes sp. isolate Gv07 small subunit ribosomal RNA gene, p	1013	1013	100%	0.0	97.98%	MF101395.1			cubensis
✓	Leucophaea lactinea strain DMC811 18S ribosomal RNA gene, p	1013	1013	100%	0.0	97.95%	KC589130.1			flavida
✓	Uncultured fungus clone CM4510 18S ribosomal RNA gene, par	1013	1013	100%	0.0	97.96%	KF800921.1			cubensis
✓	Trametes cubensis voucher T.V93_213ap 18S ribosomal RNA g	1013	1013	100%	0.0	97.96%	JN164922.1			cubensis
✓	Trametes cubensis voucher CR96 18S ribosomal RNA gene, pa	1013	1013	100%	0.0	97.96%	JN164922.1			cubensis
✓	Leucophaea lactinea strain DMC246 18S ribosomal RNA gene	1011	1011	100%	0.0	97.79%	KC589127.1			cubensis
✓	Uncultured fungus clone L042893-122-063-401 internal transcri	1011	1011	99%	0.0	99.11%	GQ999226.1			menziesii
✓	Gestrichum candidum strain C234 internal transcribed spacer 1,	1009	1009	95%	0.0	99.48%	KU377517.1			menziesii
✓	Uncultured fungus clone rcw_134 18S ribosomal RNA gene, par	1009	1009	99%	0.0	97.96%	KM104070.1			cubensis
✓	Polyporeales sp. 2 TT-2017 small subunit ribosomal RNA gene, p	1007	1007	100%	0.0	97.79%	MF621873.1			cubensis
✓	Leucophaea flavida strain DMC913 18S ribosomal RNA gene, g	1007	1007	100%	0.0	97.79%	KC589132.1			flavida
✓	Leucophaea lactinea strain DMC377 18S ribosomal RNA gene	1007	1007	100%	0.0	97.79%	KC589128.1			cubensis



Comparison of BLAST-based (pairwise alignment) vs. tree-based (multiple alignment) identification of a target fungal ITS sequence (DB42771, Vietnam; see Lücking et al. [2020](#)).

(A) A simple sequence alignment

```
AGCAATGGCCAGACAATAATG
AGCTATGGACAGACATTAATG
***  ****  *****  ****
```

(B) A more difficult sequence alignment

```
GACGACCATAGACCAGCATAG
GACTACCATAGA - CTGCAAAG
***  ****  *  ***  **
```

```
GACGACCATAGACCAGCATAG
GACTACCATAGACT - GCAAAG
***  ****  ***  **
```

Two possible positions
for the indel